

USP24 Rabbit mAb

Catalog No: #56937

Package Size: #56937-1 50ul #56937-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	USP24 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	USP24 Antibody detects endogenous levels of total USP24
Immunogen Description	A synthesized peptide derived from human USP24
Other Names	Usp24;
Accession No.	Uniprot:Q9UPU5
Uniprot	Q9UPU5
Calculated MW	294kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of USP24 expression in HeLa cell lysate.

Product Description

Involved in the ubiquitin-dependent proteolytic pathway in conjunction with the 26S proteasome.

Note: This product is for in vitro research use only