

IGF2BP2 Rabbit mAb

Catalog No: #56940

Package Size: #56940-1 50ul #56940-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	IGF2BP2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	IGF2BP2 Antibody detects endogenous levels of total IGF2BP2
Immunogen Description	A synthesized peptide derived from human IGF2BP2
Other Names	Igf2bp2; IMP2; p62; VICKZ2;
Accession No.	Uniprot:Q9Y6M1
Uniprot	Q9Y6M1
Calculated MW	66kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200

Images



Western blot analysis of IGF2BP2 expression in 293T cell lysate.

Product Description

Binds to the 5'-UTR of the insulin-like growth factor 2 (IGF2) mRNAs. Binding is isoform-specific. May regulate translation of target mRNAs.

Note: This product is for in vitro research use only