

UPP1 Rabbit mAb

Catalog No: #56941

Package Size: #56941-1 50ul #56941-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	UPP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	UPP1 Antibody detects endogenous levels of total UPP1
Immunogen Description	A synthesized peptide derived from human UPP1
Other Names	UDRPASE; UPase 1; UPASE; UPP; Upp1; UrdPase 1;
Accession No.	Uniprot:Q16831
Uniprot	Q16831
Calculated MW	28kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of UPP1 expression in (1) HeLa cell lysate; (2) 293 cell lysate.

Product Description

Catalyzes the reversible phosphorylytic cleavage of uridine and deoxyuridine to uracil and ribose- or deoxyribose-1-phosphate. The produced molecules are then utilized as carbon and energy sources or in the rescue of pyrimidine bases for nucleotide synthesis.

Note: This product is for in vitro research use only