

## ACAA2 Rabbit mAb

Catalog No: #56943

Package Size: #56943-1 50ul #56943-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	ACAA2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	ACAA2 Antibody detects endogenous levels of total ACAA2
Immunogen Description	A synthesized peptide derived from human ACAA2
Other Names	Acaa2; DSAEC;
Accession No.	Uniprot:P42765
Uniprot	P42765
Calculated MW	45kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details


WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

## Images



Western blot analysis of ACAA2 expression in HeLa cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



## Product Description

Abolishes BNIP3-mediated apoptosis and mitochondrial damage.

Note: This product is for in vitro research use only