

## UBE4A Rabbit mAb

Catalog No: #56944

Package Size: #56944-1 50ul #56944-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	UBE4A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human
Specificity	UBE4A Antibody detects endogenous levels of total UBE4A
Immunogen Description	A synthesized peptide derived from human UBE4A
Other Names	UBE4A; UBOX2; UFD2; UFD2b;
Accession No.	Uniprot:Q14139
Uniprot	Q14139
Calculated MW	123kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200FC:1:50

## Images



Western blot analysis of UBE4A expression in 293 cell lysate.

## Product Description

Binds to the ubiquitin moieties of preformed conjugates and catalyzes ubiquitin chain assembly in conjunction with E1, E2, and E3.

---

Note: This product is for in vitro research use only