

HDGF Rabbit mAb

Catalog No: #56947

Package Size: #56947-1 50ul #56947-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	HDGF Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse
Specificity	HDGF Antibody Antibody detects endogenous levels of total HDGF
Immunogen Description	A synthesized peptide derived from human HDGF
Other Names	HDGF; HMG1L2;
Accession No.	Uniprot:P51858
Uniprot	P51858
Calculated MW	37kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200FC:1:50

Images



Western blot analysis of HDGF expression in A549 cell lysate.

Product Description

Heparin-binding protein, with mitogenic activity for fibroblasts. Acts as a transcriptional repressor.

Note: This product is for in vitro research use only