

USP39 Rabbit mAb

Catalog No: #56963

Package Size: #56963-1 50ul #56963-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	USP39 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	USP39 Antibody detects endogenous levels of total USP39
Immunogen Description	A synthesized peptide derived from human USP39
Other Names	Inactive ubiquitin specific peptidase 39; SAD1; snRNP ASSEMBLY DEFECTIVE 1; SNRNP65; Ubiquitin specific peptidase 39; USP39;
Accession No.	Uniprot:Q53GS9
Uniprot	Q53GS9
Calculated MW	65kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

Images



Western blot analysis of USP39 expression in HeLa cell lysate.

Product Description

May play a role in mRNA splicing. It is unsure if the protein really exhibits hydrolase activity. Could be a competitor of ubiquitin C-terminal hydrolases (UCHs).

Note: This product is for in vitro research use only