

GLU2B Rabbit mAb

Catalog No: #56966

Package Size: #56966-1 50ul #56966-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GLU2B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	GLU2B Antibody detects endogenous levels of total GLU2B
Immunogen Description	A synthesized peptide derived from human GLU2B
Other Names	G19P1; PCLD; PKCSH; PLD1; PRKCSH;
Accession No.	Uniprot:P14314
Uniprot	P14314
Calculated MW	80kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

Images



Western blot analysis of GLU2B expression in A431 cell lysate.

Product Description

Regulatory subunit of glucosidase II.

Note: This product is for in vitro research use only