

CTNNBIP1 Rabbit mAb

Catalog No: #56967

Package Size: #56967-1 50ul #56967-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CTNNBIP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP FC
Species Reactivity	Human
Specificity	CTNNBIP1 Antibody detects endogenous levels of total CTNNBIP1
Immunogen Description	A synthesized peptide derived from human CTNNBIP1
Other Names	CTNNBIP1; ICAT;
Accession No.	Uniprot:Q9NSA3
Uniprot	Q9NSA3
Calculated MW	12kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IP:1:50FC:1:50

Images



Western blot analysis of CTNNBIP1 expression in 293T cell lysate.

Product Description

Prevents the interaction between CTNNB1 and TCF family members, and acts as negative regulator of the Wnt signaling pathway.

Note: This product is for in vitro research use only