

## COASY Rabbit mAb

Catalog No: #56968

Package Size: #56968-1 50ul #56968-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	COASY Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human
Specificity	COASY Antibody detects endogenous levels of total COASY
Immunogen Description	A synthesized peptide derived from human COASY
Other Names	CoA kinase; COASY; DPCK; DPCOAK; NBP; POV2; PPAT; UKR1;
Accession No.	Uniprot:Q13057
Uniprot	Q13057
Calculated MW	62kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details


WB:1:500~1:2000FC:1:50

## Images



Western blot analysis of COASY expression in HeLa cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



## Product Description

Bifunctional enzyme that catalyzes the fourth and fifth sequential steps of CoA biosynthetic pathway.

Note: This product is for in vitro research use only