

## TUB1 Rabbit mAb

Catalog No: #56972

Package Size: #56972-1 50ul #56972-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TUB1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	TUB1 Antibody detects endogenous levels of total TUB1
Immunogen Description	A synthesized peptide derived from human TUB1
Other Names	rd5; TUB 1; TUB;
Accession No.	Uniprot:P50607
Uniprot	P50607
Calculated MW	60kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200

## Images



Western blot analysis of TUB1 expression in Neuro-2a cell lysate.

## Product Description

Functions in signal transduction from heterotrimeric G protein-coupled receptors. Binds to membranes containing phosphatidylinositol 4,5-bisphosphate.

---

Note: This product is for in vitro research use only