

HGD Rabbit mAb

Catalog No: #56977

Package Size: #56977-1 50ul #56977-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | HGD Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB IHC ICC/IF |
| Species Reactivity | Human |
| Specificity | HGD Antibody detects endogenous levels of total HGD |
| Immunogen Description | A synthesized peptide derived from human HGD |
| Other Names | AKU; hgd; HGO; Homogentisate 1 2 dioxygenase; Homogentisate 1; Homogentisate oxidase; |
| Accession No. | Uniprot:Q93099 |
| Uniprot | Q93099 |
| Calculated MW | 48kDa |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle. |

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

Images



Western blot analysis of HGD expression in HepG2 cell lysate.

Product Description

Involved in the subpathway that synthesizes acetoacetate and fumarate from L-phenylalanine.

Note: This product is for in vitro research use only