

## ZNF259 Rabbit mAb

Catalog No: #56979

Package Size: #56979-1 50ul #56979-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	ZNF259 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	ZNF259 Antibody detects endogenous levels of total ZNF259
Immunogen Description	A synthesized peptide derived from human ZNF259
Other Names	ZNF259; ZPR1;
Accession No.	Uniprot:O75312
Uniprot	O75312
Calculated MW	51kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000 IHC: 1:50~1:200

## Images



Western blot analysis of ZNF259 expression in MCF7 cell lysate.

## Product Description

May be a signaling molecule that communicates mitogenic signals from the cytoplasm to the nucleus.

---

Note: This product is for in vitro research use only