

PDCD7 Rabbit mAb

Catalog No: #56981

Package Size: #56981-1 50ul #56981-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PDCD7 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	PDCD7 Antibody detects endogenous levels of total PDCD7
Immunogen Description	A synthesized peptide derived from human PDCD7
Other Names	ES18; HES18; Pdc7;
Accession No.	Uniprot:Q8N8D1
Uniprot	Q8N8D1
Calculated MW	55kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 IHC: 1:50~1:200 ICC/IF: 1:50~1:200

Images



Western blot analysis of PDCD7 expression in Jurkat cell lysate.

Product Description

Promotes apoptosis when overexpressed.

Note: This product is for in vitro research use only