

DNPH1 Rabbit mAb

Catalog No: #56988

Package Size: #56988-1 50ul #56988-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DNPH1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	DNPH1 Antibody detects endogenous levels of total DNPH1
Immunogen Description	A synthesized peptide derived from human DNPH1
Other Names	DNPH1; Rcl;
Accession No.	Uniprot:O43598
Uniprot	O43598
Calculated MW	21kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200

Images



Western blot analysis of DNPH1 expression in HeLa cell lysate.

Product Description

Catalyzes the cleavage of the N-glycosidic bond of deoxyribonucleoside 5'-monophosphates to yield deoxyribose 5-phosphate and a purine or pyrimidine base. Deoxyribonucleoside 5'-monophosphates containing purine bases are preferred to those containing pyrimidine bases.

Note: This product is for in vitro research use only