

Adenylate kinase 4 Rabbit mAb

Catalog No: #56989

Package Size: #56989-1 50ul #56989-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Adenylate kinase 4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	Adenylate kinase 4 Antibody detects endogenous levels of total Adenylate kinase 4
Immunogen Description	A synthesized peptide derived from human Adenylate kinase 4
Other Names	AK3; AK3L1; AK3L2; AK4;
Accession No.	Uniprot:P27144
Uniprot	P27144
Calculated MW	25kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500~1:2000 IHC: 1:50~1:200 ICC/IF: 1:50~1:200

Images



Western blot analysis of Adenylate kinase 4 expression in HepG2 cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Product Description

Catalyzes the reversible transfer of the terminal phosphate group between ATP and AMP. May also be active with GTP.

Note: This product is for in vitro research use only