

IL15RA Rabbit mAb

Catalog No: #56990



Package Size: #56990-1 50ul #56990-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description


|                       |                                                                                                                                                                                                                          |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | IL15RA Rabbit mAb                                                                                                                                                                                                        |
| Host Species          | Rabbit                                                                                                                                                                                                                   |
| Clonality             | Monoclonal                                                                                                                                                                                                               |
| Clone No.             | 72B6                                                                                                                                                                                                                     |
| Isotype               | Rabbit IgG                                                                                                                                                                                                               |
| Purification          | Affinity-chromatography                                                                                                                                                                                                  |
| Applications          | WB ICC/IF                                                                                                                                                                                                                |
| Species Reactivity    | Human                                                                                                                                                                                                                    |
| Specificity           | IL15RA Antibody detects endogenous levels of total IL15RA                                                                                                                                                                |
| Immunogen Description | A synthesized peptide derived from human IL15RA                                                                                                                                                                          |
| Target Name           | IL15RA                                                                                                                                                                                                                   |
| Other Names           | CD215; IL 15R alpha; IL-15 receptor subunit alpha; IL-15R-alpha; IL-15RA; IL15ra; Interleukin 15 receptor alpha; sIL-15 receptor subunit alpha; sIL-15R-alpha; sIL-15RA ; Soluble interleukin 15 receptor subunit alpha; |
| Accession No.         | Uniprot:Q13261                                                                                                                                                                                                           |
| Uniprot               | Q13261                                                                                                                                                                                                                   |
| Calculated MW         | 38kDa                                                                                                                                                                                                                    |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.                                                                                                                       |
| Storage               | Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.                                                                                                                                         |

Application Details

WB:1:500~1:2000

ICC/IF:1:50~1:200

Images



All Lanes: IL15RA Antibody at 1/ 500 dilution.Lane1:Raji cell lysate.Proteins at 20ug per laneSecondary: Goat Anti-Rabbit IgG(HRP) at 1/20000 dilutionPredicted band size : 38kDaObserved band size: 38kDa

## Product Description

---

Receptor for interleukin-15. Expression of different isoforms may alter or interfere with signal transduction. Isoform 5, isoform 6, isoform 7 and isoform 8 do not bind IL15. Signal transduction involves STAT3, STAT5, STAT6, JAK2 (By similarity) and SYK.

---

Note: This product is for in vitro research use only