

NCTR2 Rabbit mAb

Catalog No: #56994

Package Size: #56994-1 50ul #56994-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NCTR2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human
Specificity	NCTR2 Antibody detects endogenous levels of total NCTR2
Immunogen Description	A synthesized peptide derived from human NCTR2
Other Names	CD336; LY95; NCR2; NKp44;
Accession No.	Uniprot:O95944
Uniprot	O95944
Calculated MW	44kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 ICC/IF: 1:50~1:200

Images



Western blot analysis of NCTR2 expression in K562 cell lysate.

Product Description

Cytotoxicity-activating receptor that may contribute to the increased efficiency of activated natural killer (NK) cells to mediate tumor cell lysis.

Note: This product is for in vitro research use only