

## NDRG3 Rabbit mAb

Catalog No: #56997

Package Size: #56997-1 50ul #56997-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	NDRG3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	NDRG3 Antibody detects endogenous levels of total NDRG3
Immunogen Description	A synthesized peptide derived from human NDRG3
Other Names	NDRG3; Protein Ndr3; Protein NDRG3;
Accession No.	Uniprot:Q9UGV2
Uniprot	Q9UGV2
Calculated MW	45kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000 IHC: 1:50~1:200

## Images



Western blot analysis of NDRG3 expression in Raji cell lysate.

## Product Description

Belongs to the NDRG family.

---

Note: This product is for in vitro research use only