

HOOK2 Rabbit mAb

Catalog No: #56999

Package Size: #56999-1 50ul #56999-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	HOOK2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	HOOK2 Antibody detects endogenous levels of total HOOK2
Immunogen Description	A synthesized peptide derived from human HOOK2
Other Names	hHK2; hhook2; HK2; HOOK 2;
Accession No.	Uniprot:Q96ED9
Uniprot	Q96ED9
Calculated MW	85kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 IHC: 1:50~1:200

Images



Western blot analysis of HOOK2 expression in HeLa cell lysate.

Product Description

HOOK2 is a member of the Hook protein family. Hook proteins are cytosolic coiled-coil proteins that contain conserved N-terminal domains, which attach to microtubules, and more divergent C-terminal domains, which mediate binding to organelles.

Note: This product is for in vitro research use only