

## DPP3 Rabbit mAb

Catalog No: #57001

Package Size: #57001-1 50ul #57001-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	DPP3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human
Specificity	DPP3 Antibody detects endogenous levels of total DPP3
Immunogen Description	A synthesized peptide derived from human DPP3
Other Names	Dipeptidyl arylamidase III; Dipeptidyl peptidase 3; Dipeptidyl peptidase III; dpp3; DPPIII;
Accession No.	Uniprot:Q9NY33
Uniprot	Q9NY33
Calculated MW	83kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000FC:1:50

## Images



Western blot analysis of DPP3 expression in 293T cell lysate.

## Product Description

Cleaves Arg-Arg-beta-naphthylamide.

---

Note: This product is for in vitro research use only