

AGL Rabbit mAb

Catalog No: #57002

Package Size: #57002-1 50ul #57002-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	AGL Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	AGL Antibody detects endogenous levels of total AGL
Immunogen Description	A synthesized peptide derived from human AGL
Other Names	AGL; C77197; GDE; Glycogen debrancher;
Accession No.	Uniprot:P35573
Uniprot	P35573
Calculated MW	175kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500~1:2000

IHC: 1:50~1:200

ICC/IF: 1:50~1:200


FC: 1:50

Images



Western blot analysis of AGL expression in K562 cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Product Description

Multifunctional enzyme acting as 1,4-alpha-D-glucan:1,4-alpha-D-glucan 4-alpha-D-glycosyltransferase and amylo-1,6-glucosidase in glycogen degradation.

Note: This product is for in vitro research use only