

## NOP10 Rabbit mAb

Catalog No: #57004

Package Size: #57004-1 50ul #57004-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	NOP10 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	NOP10 Antibody detects endogenous levels of total NOP10
Immunogen Description	A synthesized peptide derived from human NOP10
Other Names	DKCB1; NOLA3; NOP10; NOP10P;
Accession No.	Uniprot:Q9NPE3
Uniprot	Q9NPE3
Calculated MW	10kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

IHC: 1:50~1:200

## Images



Western blot analysis of NOP10 expression in Jurkat cell lysate.

## Product Description

Required for ribosome biogenesis and telomere maintenance. Part of the H/ACA small nucleolar ribonucleoprotein (H/ACA snoRNP) complex, which catalyzes pseudouridylation of rRNA.

---

Note: This product is for in vitro research use only