

NPY5R Rabbit mAb

Catalog No: #57005

Package Size: #57005-1 50ul #57005-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NPY5R Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human
Specificity	NPY5R Antibody detects endogenous levels of total NPY5R
Immunogen Description	A synthesized peptide derived from human NPY5R
Other Names	NPY 5R; NPY R5; NPY Y5; NPY Y5 receptor; NPY5R; NPYR5; NPYY5;Y5 receptor;
Accession No.	Uniprot:Q15761
Uniprot	Q15761
Calculated MW	51kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200FC:1:50

Images



Western blot analysis of NPY5R expression in Ramos cell lysate.

Product Description

Receptor for neuropeptide Y and peptide YY. The activity of this receptor is mediated by G proteins that inhibit adenylate cyclase activity. Seems to be associated with food intake. Could be involved in feeding disorders.

Note: This product is for in vitro research use only