

## DCAF7 Rabbit mAb

Catalog No: #57006

Package Size: #57006-1 50ul #57006-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	DCAF7 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC FC
Species Reactivity	Human
Specificity	DCAF7 Antibody detects endogenous levels of total DCAF7
Immunogen Description	A synthesized peptide derived from human DCAF7
Other Names	AN11; Dcaf7; HAN11; WDR68;
Accession No.	Uniprot:P61962
Uniprot	P61962
Calculated MW	39kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200FC:1:50

## Images



Western blot analysis of DCAF7 expression in HeLa cell lysate.

## Product Description

Involved in craniofacial development. Acts upstream of the EDN1 pathway and is required for formation of the upper jaw equivalent, the palatoquadrate.

---

Note: This product is for in vitro research use only