

TAF15 Rabbit mAb

Catalog No: #57011

Package Size: #57011-1 50ul #57011-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TAF15 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	TAF15 Antibody detects endogenous levels of total TAF15
Immunogen Description	A synthesized peptide derived from human TAF15
Other Names	Npl3; RBP56; RNA binding protein 56; TAF; TAF15; TAF2N; TAFII68;
Accession No.	Uniprot:Q92804
Uniprot	Q92804
Calculated MW	75kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200FC:1:50

Images



Western blot analysis of TAF15 expression in HeLa cell lysate.

Product Description

RNA and ssDNA-binding protein that may play specific roles during transcription initiation at distinct promoters. Can enter the preinitiation complex together with the RNA polymerase II (Pol II).

Note: This product is for in vitro research use only