

CNAP1 Rabbit mAb

Catalog No: #57017

Package Size: #57017-1 50ul #57017-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CNAP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	CNAP1 Antibody detects endogenous levels of total CNAP1
Immunogen Description	A synthesized peptide derived from human CNAP1
Other Names	CAPD2; CEP2; CEP250; CNAP1; CND1; Condensin; hCAPD2; NCAPD2;
Accession No.	Uniprot:Q15021
Uniprot	Q15021
Calculated MW	150kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of CNAP1 expression in HeLa cell lysate.

Product Description

CNAP1/hCAP-D2 is a regulatory subunit of the condensin complex, a complex required for conversion of interphase chromatin into mitotic-like condensed chromosomes.

Note: This product is for in vitro research use only