

Thromboxane A2 receptor Rabbit mAb

Catalog No: #57022

Package Size: #57022-1 50ul #57022-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Thromboxane A2 receptor Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Thromboxane A2 receptor Antibody detects endogenous levels of total Thromboxane A2 receptor
Immunogen Description	A synthesized peptide derived from human Thromboxane A2 receptor
Other Names	BDPLT13; Tbx2r; TP receptor; TXA2R;
Accession No.	Uniprot:P21731
Uniprot	P21731
Calculated MW	50kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of Thromboxane A2 receptor expression in human platelet cell lysate.

Product Description

Receptor for thromboxane A2 (TXA2), a potent stimulator of platelet aggregation.

Note: This product is for in vitro research use only