

## RPL13 Rabbit mAb

Catalog No: #57023

Package Size: #57023-1 50ul #57023-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	RPL13 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse
Specificity	RPL13 Antibody detects endogenous levels of total RPL13
Immunogen Description	A synthesized peptide derived from human RPL13
Other Names	BBC1; L13; rpl13;
Accession No.	Uniprot:P26373
Uniprot	P26373
Calculated MW	26kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

## Images



Western blot analysis of RPL13 expression in MCF7 cell lysate.

## Product Description

Belongs to the ribosomal protein L13e family.

---

Note: This product is for in vitro research use only