

## TRPV5 Rabbit mAb

Catalog No: #57024

Package Size: #57024-1 50ul #57024-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TRPV5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	TRPV5 Antibody detects endogenous levels of total TRPV5
Immunogen Description	A synthesized peptide derived from human TRPV5
Other Names	CAT2; ECaC; ECaC1; OTRPC3; TrpV5;
Accession No.	Uniprot:Q9NQA5
Uniprot	Q9NQA5
Calculated MW	90kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000 IHC: 1:50~1:200

## Images



Western blot analysis of TRPV5 expression in HepG2 cell lysate.

## Product Description

Constitutively active calcium selective cation channel thought to be involved in Ca(2+) reabsorption in kidney and intestine.

---

Note: This product is for in vitro research use only