

PKN2 Rabbit mAb

Catalog No: #57026

Package Size: #57026-1 50ul #57026-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PKN2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	PKN2 Antibody detects endogenous levels of total PKN2
Immunogen Description	A synthesized peptide derived from human PKN2
Other Names	PAK2; PKN gamma; PKN2; PRK2; PRKCL2;
Accession No.	Uniprot:Q16513
Uniprot	Q16513
Calculated MW	140kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 ICC/IF: 1:50~1:200

Images



Western blot analysis of PKN2 expression in HeLa cell lysate.

Product Description

Exhibits a preference for highly basic protein substrates.

Note: This product is for in vitro research use only