

HARS Rabbit mAb

Catalog No: #57028

Package Size: #57028-1 50ul #57028-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | HARS Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB IHC |
| Species Reactivity | Human |
| Specificity | HARS Antibody detects endogenous levels of total HARS |
| Immunogen Description | A synthesized peptide derived from human HARS |
| Other Names | HARS; HisRS; histidine transylase; HRS; USH3B; |
| Accession No. | Uniprot:P12081 |
| Uniprot | P12081 |
| Calculated MW | 50kDa |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle. |

Application Details

WB: 1:500~1:2000 IHC: 1:50~1:200

Images



Western blot analysis of HARS expression in 293T cell lysate.

Product Description

Plays a role in axon guidance.

Note: This product is for in vitro research use only