

RPL26L1 Rabbit mAb

Catalog No: #57032

Package Size: #57032-1 50ul #57032-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	RPL26L1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	RPL26L1 Antibody detects endogenous levels of total RPL26L1
Immunogen Description	A synthesized peptide derived from human RPL26L1
Other Names	RPL26P1;
Accession No.	Uniprot:Q9UNX3
Uniprot	Q9UNX3
Calculated MW	17kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

Images



Western blot analysis of RPL26L1 expression in HepG2 cell lysate.

Product Description

Belongs to the ribosomal protein L24P family.

Note: This product is for in vitro research use only