

PCBD1 Rabbit mAb

Catalog No: #57035

Package Size: #57035-1 50ul #57035-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PCBD1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	PCBD1 Antibody detects endogenous levels of total PCBD1
Immunogen Description	A synthesized peptide derived from human PCBD1
Other Names	DCoH; PCBD; PCBD1; PCD; PHS;
Accession No.	Uniprot:P61457
Uniprot	P61457
Calculated MW	12kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200

Images



Western blot analysis of PCBD1 expression in Caco-2 cell lysate.

Product Description

Involved in tetrahydrobiopterin biosynthesis. Seems to both prevent the formation of 7-pterins and accelerate the formation of quinonoid-BH2. Coactivator for HNF1A-dependent transcription.

Note: This product is for in vitro research use only