

COPS4 Rabbit mAb

Catalog No: #57037

Package Size: #57037-1 50ul #57037-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	COPS4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	COPS4 Antibody detects endogenous levels of total COPS4
Immunogen Description	A synthesized peptide derived from human COPS4
Other Names	CH4; COPS4; CSN4; DCH4; SGN4;
Accession No.	Uniprot:Q9BT78
Uniprot	Q9BT78
Calculated MW	50kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of COPS4 expression in SH-SY5Y cell lysate.

Product Description

Component of the COP9 signalosome complex (CSN), a complex involved in various cellular and developmental processes.

Note: This product is for in vitro research use only