

PNMT Rabbit mAb

Catalog No: #57040

Package Size: #57040-1 50ul #57040-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PNMT Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human Mouse
Specificity	PNMT Antibody detects endogenous levels of total PNMT
Immunogen Description	A synthesized peptide derived from human PNMT
Other Names	PENT; PNMT; PNMTase;
Accession No.	Uniprot:P11086
Uniprot	P11086
Calculated MW	30kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 ICC/IF: 1:50~1:200

Images



Western blot analysis of PNMT expression in K562 cell lysate.

Product Description

Converts noradrenaline to adrenaline.

Note: This product is for in vitro research use only