

P2Y12 Rabbit mAb

Catalog No: #57043

Package Size: #57043-1 50ul #57043-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	P2Y12 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	P2Y12 Antibody detects endogenous levels of total P2Y12
Immunogen Description	A synthesized peptide derived from human P2Y12
Other Names	ADP glucose receptor; ADPGR; BDPLT8; HORK3; P2RY12; P2Y purinoceptor 12; P2Y12; P2Y12 platelet ADP receptor; Platelet ADP receptor; Purinergic receptor P2RY12; Purinergic receptor P2Y12; SP1999;
Accession No.	Uniprot:Q9H244
Uniprot	Q9H244
Calculated MW	50kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

IP:1:50

Images



Western blot analysis of P2Y12 expression in RAW 264.7 cell lysate.

Product Description

Receptor for ADP and ATP coupled to G-proteins that inhibit the adenylyl cyclase second messenger system. Not activated by UDP and UTP. Involved in platelets aggregation.

Note: This product is for in vitro research use only