

HNF-4-alpha Rabbit mAb

Catalog No: #57044

Package Size: #57044-1 50ul #57044-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | HNF-4-alpha Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB IHC IP CHIP |
| Species Reactivity | Human Mouse Rat |
| Specificity | HNF-4-alpha Antibody detects endogenous levels of total HNF-4-alpha |
| Immunogen Description | A synthesized peptide derived from human HNF-4-alpha |
| Other Names | FRTS4; Hepatic nuclear factor 4 alpha; Hepatocyte nuclear factor 4 alpha; Hepatocyte nuclear factor 4; HNF 4 alpha; HNF 4; HNF4; HNF4A; HNF4a7; HNF4a8; HNF4a9; Hnf4alpha; HNF4alpha10/11/12; MODY 1; MODY; MODY1; NR2A1; NR2A21; TCF; TCF14; Tcf4; Transcription factor 14; Transcription factor HNF4; |
| Accession No. | Uniprot:P41235 |
| Uniprot | P41235 |
| Calculated MW | 53kDa |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle. |

Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:50

Images



Western blot analysis of HNF-4-alpha expression in HepG2 cell lysate.

Product Description

Transcriptionally controlled transcription factor. Binds to DNA sites required for the transcription of alpha 1-antitrypsin, apolipoprotein CIII, transthyretin genes and HNF1-alpha. May be essential for development of the liver, kidney and intestine.

Note: This product is for in vitro research use only