

GFPT1 Rabbit mAb

Catalog No: #57052

Package Size: #57052-1 50ul #57052-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GFPT1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Mouse Rat
Specificity	GFPT1 Antibody detects endogenous levels of total GFPT1
Immunogen Description	A synthesized peptide derived from human GFPT1
Other Names	CMS12; CMSTA1; GFA; GFAT; GFAT1; GFAT1m; GFPT; Gfpt1; GFPT1L; Glucosamine--fructose-6-phosphate aminotransferase [isomerizing] 1; Hexosephosphate aminotransferase 1;
Accession No.	Uniprot:Q06210
Uniprot	Q06210
Calculated MW	79kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:50

Images



Western blot analysis of GFPT1 expression in MCF7 cell lysate.

Product Description

Controls the flux of glucose into the hexosamine pathway. Most likely involved in regulating the availability of precursors for N- and O-linked glycosylation of proteins.

Note: This product is for in vitro research use only