

## FRAT1 Rabbit mAb

Catalog No: #57061

Package Size: #57061-1 50ul #57061-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	FRAT1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	FRAT1 Antibody detects endogenous levels of total FRAT1
Immunogen Description	A synthesized peptide derived from human FRAT1
Other Names	FRAT 1;
Accession No.	Uniprot:Q92837
Uniprot	Q92837
Calculated MW	29kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000 ICC/IF: 1:50~1:200

## Images



Western blot analysis of FRAT1 expression in HeLa cell lysate.

## Product Description

Positively regulates the Wnt signaling pathway by stabilizing beta-catenin through the association with GSK-3. May play a role in tumor progression and collaborate with PIM1 and MYC in lymphomagenesis.

---

Note: This product is for in vitro research use only