

## MVK Rabbit mAb

Catalog No: #57072

Package Size: #57072-1 50ul #57072-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	MVK Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human
Specificity	MVK Antibody detects endogenous levels of total MVK
Immunogen Description	A synthesized peptide derived from human MVK
Other Names	LRBP; Mevalonate kinase 1; mvk; MVLK; POROK3;
Accession No.	Uniprot:Q03426
Uniprot	Q03426
Calculated MW	42kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000FC:1:50

## Images



Western blot analysis of MVK expression in HepG2 cell lysate.

## Product Description

May be a regulatory site in cholesterol biosynthetic pathway.

---

Note: This product is for in vitro research use only