

PPCS Rabbit mAb

Catalog No: #57077

Package Size: #57077-1 50ul #57077-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PPCS Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	PPCS Antibody detects endogenous levels of total PPCS
Immunogen Description	A synthesized peptide derived from human PPCS
Other Names	COAB; PPC synthetase; Ppcs;
Accession No.	Uniprot:Q9HAB8
Uniprot	Q9HAB8
Calculated MW	34kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 IHC: 1:50~1:200 ICC/IF: 1:50~1:200

Images



Western blot analysis of PPCS expression in SH-SY5Y cell lysate.

Product Description

Catalyzes the first step in the biosynthesis of coenzyme A from vitamin B5, where cysteine is conjugated to 4'-phosphopantothenate to form 4-phosphopantothenoylecysteine.

Note: This product is for in vitro research use only