

NUP50 Rabbit mAb

Catalog No: #57078

Package Size: #57078-1 50ul #57078-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NUP50 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	NUP50 Antibody detects endogenous levels of total NUP50
Immunogen Description	A synthesized peptide derived from human NUP50
Other Names	NPAP60; NPAP60L; Nucleoporin 50kDa; Nup50;
Accession No.	Uniprot:Q9UKX7
Uniprot	Q9UKX7
Calculated MW	55kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 ICC/IF: 1:50~1:200

Images



Western blot analysis of NUP50 expression in HeLa cell lysate.

Product Description

Component of the nuclear pore complex that has a direct role in nuclear protein import. Actively displaces NLSs from importin-alpha, and facilitates disassembly of the importin-alpha:beta-cargo complex and importin recycling.

Note: This product is for in vitro research use only