

## DTYMK Rabbit mAb

Catalog No: #57081

Package Size: #57081-1 50ul #57081-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	DTYMK Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	DTYMK Antibody detects endogenous levels of total DTYMK
Immunogen Description	A synthesized peptide derived from human DTYMK
Other Names	CDC8; Dtymk; TMPK; TYMK;
Accession No.	Uniprot:P23919
Uniprot	P23919
Calculated MW	23kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200FC:1:50

## Images



Western blot analysis of DTYMK expression in HeLa cell lysate.

## Product Description

Catalyzes the conversion of dTMP to dTDP.

---

Note: This product is for in vitro research use only