

TPMT Rabbit mAb

Catalog No: #57082

Package Size: #57082-1 50ul #57082-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TPMT Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC FC
Species Reactivity	Human Mouse Rat
Specificity	TPMT Antibody detects endogenous levels of total TPMT
Immunogen Description	A synthesized peptide derived from human TPMT
Other Names	Thiopurine methyltransferase; TPMT;
Accession No.	Uniprot:P51580
Uniprot	P51580
Calculated MW	30kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 IHC: 1:50~1:200 FC: 1:50

Images



Western blot analysis of TPMT expression in Jurkat cell lysate.

Product Description

Catalyzes the S-methylation of thiopurine drugs such as 6-mercaptopurine.

Note: This product is for in vitro research use only