

MAS1L Rabbit mAb

Catalog No: #57084

Package Size: #57084-1 50ul #57084-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MAS1L Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human
Specificity	MAS1L Antibody detects endogenous levels of total MAS1L
Immunogen Description	A synthesized peptide derived from human MAS1L
Other Names	MAS L; MAS R; MAS1L; MRG;
Accession No.	Uniprot:P35410
Uniprot	P35410
Calculated MW	37kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 IP: 1:50

Images



Western blot analysis of MAS1L expression in 293T cell lysate.

Product Description

Belongs to the G-protein coupled receptor 1 family. Mas subfamily.

Note: This product is for in vitro research use only