

PNK / PNKP Rabbit mAb

Catalog No: #57087

Package Size: #57087-1 50ul #57087-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PNK / PNKP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	PNK / PNKP Antibody detects endogenous levels of total PNK / PNKP
Immunogen Description	A synthesized peptide derived from human PNK / PNKP
Other Names	DEM1; EIEE10; MCSZ; PNK1; Pnkp;
Accession No.	Uniprot:Q96T60
Uniprot	Q96T60
Calculated MW	57kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

Images



Western blot analysis of PNK / PNKP expression in Jurkat cell lysate.

Product Description

Catalyzes the phosphorylation of DNA at 5'-hydroxyl termini and can dephosphorylate its 3'-phosphate termini. Plays an important function in DNA repair following ionizing radiation or oxidative damage.

Note: This product is for in vitro research use only