

RAVER2 Rabbit mAb

Catalog No: #57091

Package Size: #57091-1 50ul #57091-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	RAVER2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human
Specificity	RAVER2 Antibody detects endogenous levels of total RAVER2
Immunogen Description	A synthesized peptide derived from human RAVER2
Other Names	Protein raver 2; Raver2;
Accession No.	Uniprot:Q9HCJ3
Uniprot	Q9HCJ3
Calculated MW	74kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50

Images



Western blot analysis of RAVER2 expression in 293T cell lysate.

Product Description

May bind single-stranded nucleic acids.

Note: This product is for in vitro research use only