

NDUFA1 Rabbit mAb

Catalog No: #57092

Package Size: #57092-1 50ul #57092-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NDUFA1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	NDUFA1 Antibody detects endogenous levels of total NDUFA1
Immunogen Description	A synthesized peptide derived from human NDUFA1
Other Names	NDUFA1; ZNF183;
Accession No.	Uniprot:O15239
Uniprot	O15239
Calculated MW	10kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 IHC: 1:50~1:200

Images



Western blot analysis of NDUFA1 expression in A673 cell lysate.

Product Description

Accessory subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I), that is believed not to be involved in catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.

Note: This product is for in vitro research use only