

SFT Rabbit mAb

Catalog No: #57096

Package Size: #57096-1 50ul #57096-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SFT Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Mouse Rat
Specificity	SFT Antibody detects endogenous levels of total SFT
Immunogen Description	A synthesized peptide derived from human SFT
Other Names	SFT; UBC4/5; UbcH5; UBCH5A; Ube2d1;
Accession No.	Uniprot:P51668
Uniprot	P51668
Calculated MW	17kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:50

Images



Western blot analysis of SFT expression in Raji cell lysate.

Product Description

Accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. In vitro catalyzes 'Lys-11'-linked polyubiquitination. Involved in the selective degradation of short-lived and abnormal proteins. Down-regulated during the S-phase it is involved in progression through the cell cycle.

Note: This product is for in vitro research use only